
My COVID-19 test is expiring! What to know.

Check the box for the final expiration date for this test

How to find the control line






Add liquid sample.
Control line should appear. If not, discard test.



*Time will vary depending on the brand of the COVID-19 self-test

If positive, the test line will appear within 10-30* minutes. Positive results may appear bold or faint

or faint.



Trust your positive results!

A false positive is highly unlikely.

You can use an expired test as long as the control line is working.



For more COVID-19 testing information call the ACPHD Front Door Line at 510-577-7080 or visit covid-19.acgov.org/testing

My COVID-19 test is expiring! What to know.

Check the box for the final expiration date for this test

How to find the control line



How the test should look out of the box



Add liquid sample.
Control line should appear. If not, discard test.



*Time will vary depending on the brand of the COVID-19 self-test

If positive, the test line will appear within 10-30* minutes.
Positive results may appear bold

or faint.



Trust your positive results!

A false positive is highly unlikely.

You can use an expired test as long as the control line is working.

